

piccolo[®]

STAT Chemistry Analyzer

piccolo[®]
Panels

CLIA
Waived

Lipid Panel Plus

ALT, AST, CHOL, CHOL/HDL*
GLU, HDL, LDL*, TRIG, VLDL*

Liver Panel Plus

ALB, ALP, ALT, AMY
AST, CGT, TBIL, TP

Electrolyte

Cl⁻, K⁺, Na⁺, tCO₂

Metlyte 8

BUN, CK, Cl⁻, CRE
GLU, K⁺, Na⁺, tCO₂

CLIA
Waived

Lipid

CHOL, CHOL/HDL*, HDL
LDL*, TRIG, VLDL*

Basic Metabolic

BUN, Ca, Cl⁻, CRE
GLU, K⁺, Na⁺, tCO₂

Renal Function

ALB, BUN, Ca, Cl⁻, CRE
GLU, K⁺, Na⁺, PHOS, tCO₂

Hepatic Function

ALB, ALP, ALT, AST
DBIL, TBIL, TP

General Chemistry 6

ALT, AST, BUN, CRE, CGT, GLU

General Chemistry 13

ALB, ALP, ALT, AMY, AST, BUN
Ca, CRE, CGT, GLU, TBIL, TP, UA

Comprehensive Metabolic

ALB, ALP, ALT, AST, BUN, Ca, Cl⁻
CRE, GLU, K⁺, Na⁺, TBIL, tCO₂, TP

Basic Metabolic Panel Plus

BUN, Ca, Cl⁻, CRE, GLU, K⁺, Na⁺
tCO₂, Mg, Lactate Dehydrogenase

*calculated

ABAXIS



piccolo[®]

STAT Chemistry Analyzer

Advanced Point-of-Care Chemistry System

Piccolo Features

As Easy as 1-2-3

1. Pipette a 100 µl sample (whole blood, serum or plasma) into the self-contained reagent disc
2. Place the disc in the Piccolo and enter basic patient information
3. Get printed results in 12 minutes.

No Lab Tech Required

As easy to run as a CD player, anyone in the office can run the Piccolo

Self Calibrating

Each Reagent Disc is barcoded; No manual calibration ever.

Only Run Controls Every 30 Days

On-Board intelligent Quality Controls meet CLIA and NCCLS requirements for daily control testing

Extensive Test Menu!

The Piccolo is the only instrument offering the most commonly requested chemistry panels

Smallest Footprint Available

Compact and portable, the Piccolo fits perfectly in labs with limited space! (H9.5" h X W6" w X D11.5", 15 lbs.)

CLIA Waived Lipid Testing!

Abaxis assists you in meeting CLIA requirements for Non-Waived testing

Lab Accurate Results

Correlation studies with Hospital and Commercial labs show >96% accuracy

Space-Age Technology

Patented design originally created for NASA and used by the U.S. Military

Practice Impact

Improves Patient Care

Rapid response to abnormal labs allows for instant diagnosis, consultation, treatment and real time adjustments to medication while the patient is in the office

Enhances Patient Satisfaction

Treating the patient in one visit increases the value of your services to your client

Improves Practice Efficiencies

Reduces staff time chasing charts, contacting patients regarding delayed results, or rescheduling the patient back for another visit

May Lower Liability Risk

Immediate results may lessen medical error or negligence liability.

Easy Practice Integration

The Piccolo's small footprint, ease of use, and quick results, fits nicely into current practice routines

Financial Impact

Increase CPT Reimbursement

The Piccolo can recover the revenue stream from the commercial lab and put it back into your practice! All Piccolo Panels are CMS Reimbursable

Enhanced E/M Billing

The time spent discussing Point-of-Care lab results with patients during the office visit enables higher level E/M Coding

No Wasted Reagents!

Since each self-contained reagent disc has a lengthy refrigerated shelf life, run the panel when you need them and worry less about expired product

No Maintenance Cost

The Piccolo's solid state design requires no pumps, tubes, or moving parts (other than the built-in centrifuge)

ABAXIS Lead Qualification Form

1. Account Name _____ Date _____

2. Contact Name _____ Title _____

3. Address _____ City _____

4. State _____ Zip _____ Phone _____ Fax: _____

5. E-mail _____

6. Distributor _____ Representative _____

7. How many chemistry panels do you request each day? _____

8. Number of Panels per day: Lipid _____ CMP _____ BMP _____
Hepatic _____ Renal _____ BUN/Creat _____ Other _____

9. Current CLIA status level (check one): Waived Non-Waived (Moderate)

Fax completed form to (510) 892-2995 to schedule a demonstration.